

AFIAS Tumormarker Control

Lot TMCOOA05

Verfall: 2019-01



PSA Level 1

Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	23%
Zielwert	0.91 ng/ml
Standardabweichung	0.07 ng/ml

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
0.70	0.77	0.84	0.91	0.98	1.05	1.12

PSA Level 2

Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert	9.46 ng/ml
Standardabweichung	0.75 ng/ml

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
7.21	7.96	8.71	9.46	10.21	10.96	11.71

--	--
