

# Liquid Diabetes Control Cliniq

## Skyla Hi

LOT 1706085/1706086A+B

Exp. 2019-02-28



### HbA1c

Level 1

Max. Toleranzbreite QUALAB	9%
Toleranzbreite Hersteller	9%
Zielwert (%)	6.9
Standardabweichung (%)	0.2

Einheit	mg/l
Datum	Messwert


-3s	-2s	-1s	Mw	+1s	+2s	+3s
6.3	6.5	6.7	6.9	7.1	7.3	7.5


### HbA1c

Level 2

Max. Toleranzbreite QUALAB	9%
Toleranzbreite Hersteller	9%
Zielwert (%)	13.1
Standardabweichung (%)	0.4

Einheit	mg/l
Datum	Messwert


-3s	-2s	-1s	Mw	+1s	+2s	+3s
11.9	12.3	12.7	13.1	13.5	13.9	14.3
