

LOT		MX 453 L		CONTROL		L		(Exp.)		2025-08-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST				60 8P - 60 18P					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	2.1				2.1					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.24				2.23					± 0.17
HB HGB	g/dl	6.2				6.2					± 0.5
	g/l	62				62					± 5
	mmol/l	3.85				3.85					± 0.31
HT HCT	%	16.6				16.5					± 2.0
	l/l	0.166				0.165					± 0.020
VGM MCV	µm ³ ; fl	74				74					± 4
TGMH MCH	pg	27.7				27.8					± 2.0
	fmol	1.72				1.73					± 0.12
CCMH MCHC	g/dl	37.4				37.6					± 3.0
	g/l	374				376					± 30
	mmol/l	23.23				23.33					± 1.86
IDR RDW	%	16.0				14.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	87				84					± 22
VPM MPV	µm ³ ; fl	7.7				7.0					± 2.0
LYMPHO	#	1.31				1.30					± 0.40
	%	62.5				62.0					± 8.0
MONO	#	0.21				0.22					± 0.19
	%	10.0				10.5					± 9.0
GRANULO	#	0.58				0.58					± 0.40
	%	27.5				27.5					± 10.0

LOT		MX 453 N		CONTROL		N		 (Exp.)		2025-08-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST				60 8P - 60 18P					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.7				7.6					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.69				4.72					± 0.26
HB HGB	g/dl	13.6				13.8					± 0.7
	g/l	136				138					± 7
	mmol/l	8.45				8.57					± 0.43
HT HCT	%	38.0				37.8					± 2.7
	l/l	0.380				0.378					± 0.027
VGM MCV	µm ³ ; fl	81				80					± 4
TGMH MCH	pg	29.0				29.2					± 2.0
	fmol	1.80				1.82					± 0.12
CCMH MCHC	g/dl	35.8				36.5					± 3.0
	g/l	358				365					± 30
	mmol/l	22.23				22.70					± 1.86
IDR RDW	%	15.5				14.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	273				263					± 40
VPM MPV	µm ³ ; fl	7.2				6.7					± 2.0
LYMPHO	#	2.46				2.47					± 0.40
	%	32.0				32.5					± 5.0
MONO	#	0.50				0.49					± 0.44
	%	6.5				6.5					± 5.9
GRANULO	#	4.74				4.64					± 0.76
	%	61.5				61.0					± 10.0

Rel: TEMP-0830 Rev.43 FRONT / RECTO 1300104144

LOT		MX 453 H		Rev 1		CONTROL		H		(Exp.)		2025-08-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST				60 8P - 60 18P							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	21.0				20.9						± 1.7	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.76				5.80						± 0.32	
HB HGB	g/dl	18.2				18.5						± 0.9	
	g/l	182				185						± 9	
	mmol/l	11.30				11.49						± 0.56	
HT HCT	%	51.3				51.6						± 3.4	
	l/l	0.513				0.516						± 0.034	
VGM MCV	µm ³ ; fl	89				89						± 4	
TGMH MCH	pg	31.6				31.9						± 2.0	
	fmol	1.96				1.98						± 0.12	
CCMH MCHC	g/dl	35.5				35.8						± 3.0	
	g/l	355				358						± 30	
	mmol/l	22.05				22.26						± 1.86	
IDR RDW	%	15.0				14.0						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	529				515						± 65	
VPM MPV	µm ³ ; fl	7.0				6.5						± 2.0	
LYMPHO	#	3.36				3.34						± 1.00	
	%	16.0				16.0						± 5.0	
MONO	#	0.84				0.84						± 0.76	
	%	4.0				4.0						± 3.6	
GRANULO	#	16.80				16.72						± 2.09	
	%	80.0				80.0						± 10.0	