

AFIAS Tumormarker Control

Lot TMC00F06

Verfall: 2019.06



PSA

Level 1

Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert	1.10 ng/ml
Standardabweichung	0.08 ng/ml

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
0.86	0.94	1.02	1.10	1.18	1.26	1.34

PSA

Level 2

Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert	9.77 ng/ml
Standardabweichung	0.78 ng/ml

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
7.43	8.21	8.99	9.77	10.55	11.33	12.11

--	--

--	--	--	--	--	--	--	--	--	--
