

AFIAS Ferritin Control

Lot FRCOOG07 Verfall: 2019.07

Ferritin Level 1

Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert	57.0 ng/ml
Standardabweichung	4.5 ng/ml

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
43.5	48.0	52.5	57.0	61.5	66.0	70.5

Ferritin Level 2

Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert	170.4 ng/ml
Standardabweichung	13.6 ng/ml

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
129.6	143.2	156.8	170.4	184.0	197.6	211.2

--	--

--	--	--	--	--	--	--	--
