

AFIAS Ferritin Control

LOT FRCOOK08 Exp. 2019.11



Ferritin

Level 1

Max. Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert (ng/ml)	63.60
Standardabweichung (ng/ml)	5.08

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
48.36	53.44	58.52	63.60	68.68	73.76	78.84

Ferritin

Level 2

Max. Toleranzbreite QUALAB	24%
Toleranzbreite Hersteller	24%
Zielwert (ng/ml)	202.2
Standardabweichung (ng/ml)	16.17

Einheit	ng/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
153.69	169.86	186.03	202.20	218.37	234.54	250.71