

AFIAS Hormone Control

LOT HOCOPH13 Exp. 2020.08



TSH

Level 1

Max. Toleranzbreite QUALAB	20%
Toleranzbreite Hersteller	20%
Zielwert (uIU/ml)	1.10
Standardabweichung (uIU/ml)	0.07

Einheit	uIU/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
0.89	0.96	1.03	1.10	1.17	1.24	1.31

TSH

Level 2

Max. Toleranzbreite QUALAB	20%
Toleranzbreite Hersteller	20%
Zielwert (uIU/ml)	7.9
Standardabweichung (uIU/ml)	0.52

Einheit	uIU/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
6.34	6.86	7.38	7.90	8.42	8.94	9.46

AFIAS Hormone Control

LOT HOCOPH13 Exp. 2020.08



T3

Level 1

Toleranzbreite QUALAB	20%
Toleranzbreite Hersteller	20%
Zielwert (nmol/L)	1.40
Standardabweichung (nmol/L)	0.09

Einheit	nmol/L
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
1.13	1.22	1.31	1.40	1.49	1.58	1.67

T3

Level 2

Toleranzbreite QUALAB	20%
Toleranzbreite Hersteller	20%
Zielwert (nmol/L)	4.10
Standardabweichung (nmol/L)	0.27

Einheit	nmol/L
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
3.29	3.56	3.83	4.10	4.37	4.64	4.91

AFIAS Hormone Control

LOT HOCOPH13 Exp. 2020.08



βhCG

Level 1

Toleranzbreite QUALAB	25%
Toleranzbreite Hersteller	25%
Zielwert (mIU/ml)	177.9
Standardabweichung (mIU/ml)	14.8

Einheit	mIU/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
133.5	148.3	163.1	177.9	192.7	207.5	222.3

βhCG

Level 2

Toleranzbreite QUALAB	25%
Toleranzbreite Hersteller	25%
Zielwert (mIU/ml)	2123.30
Standardabweichung (mIU/ml)	176.9

Einheit	mIU/ml
Datum	Messwert

-3s	-2s	-1s	Mw	+1s	+2s	+3s
1592.6	1769.5	1946.4	2123.3	2300.2	2477.1	2654.0

